

# NeuroGuard




*Cutting Edge Precision without Cutting Nerves*

*ECE 445 Senior Design Course Launch, 8/26/25*

# Our Team




Nikita Chigullapally  
CIMED M2  
*nikitc3@illinois.edu*



Prateek Dullur  
CIMED M2  
*pdullur2@illinois.edu*



Meenakshi Singhal  
CIMED M2  
*ms189@illinois.edu*



Nicolas Kelhofer  
CIMED M2  
*nck7@illinois.edu*



Kathure Kiunga  
Undergraduate  
*skiunga2@illinois.edu*

## Clinical Advisors



Suguna Pappu, MD, PhD  
*Carle Health, Neurosurgeon*



Victor Stams, MD  
*Carle Health, Plastic surgeon*



Yuan Yang, PhD  
*UIUC Bioengineering*



Michael Oelze, PhD  
*UIUC Electrical Engineering*




Arijit Banerjee, PhD  
*UIUC Electrical Engineering*

## Engineering Advisors


# 40-60%

## Nerve Injury rate with mastectomy


## Over 100,000/year in US



# An Electrical Engineering Problem, not a clinical one




# Electrocautery components



- Foot pedals
- Cautery unit
- Probe
- Grounding pad (monopolar)

# Energy Delivery by Duty Cycle



## OUTPUT CHARACTERISTICS



### Maximum Output for Monopolar and Bipolar Modes

Power readouts agree with actual power into rated load to within 20% or 5 watts, whichever is greater.

Mode	Output Power	Output Frequency	Repetition Rate	Open Circuit V <sub>peak max</sub>	Crest Factor* (Rated Load)
Cut	120 W @ 500 $\Omega$	357 kHz $\pm$ 50 kHz	N / A	1250V	2.9 $\pm$ 20%
Blend	90 W @ 800 $\Omega$	357 kHz $\pm$ 50 kHz	30 kHz $\pm$ 5 kHz	1850V	3.3 $\pm$ 20%
Coagulation	80 W @ 1000 $\Omega$	475 kHz $\pm$ 19 kHz	57 kHz $\pm$ 5 kHz	3300V	5.5 $\pm$ 20%
Fulguration	40 W @ 1000 $\Omega$	410 kHz $\pm$ 50 kHz	25 kHz $\pm$ 5 kHz	3900V	7.7 $\pm$ 20%
Bipolar	30 W @ 200 $\Omega$	520 kHz (-14 kHz, +29 kHz)	32 kHz $\pm$ 5 kHz	1200V	6.9 $\pm$ 20%



\* an indication of a waveform's ability to coagulate bleeders without a cutting effect

# An Electrical Engineering Problem, not a clinical one





# Neuroguard Components



Electrocautery waveform (~1000 V, 510 kHz)


Nerve stimulation waveform (~0.1 V, 100 Hz)








# Integration with existing tools

- We envision our device as an add on to electrocautery generator
  - Lower costs, more accessible, small learning curve
  - NeuroGuard probes which are compatible with Bovie and smoke evac systems



# Prelim Circuit Design



Current iteration will be provided to participating senior design team

# ECE 445, Fall 2025 Project

## Deliverables

- A. Simulation of circuit demonstrating ability to switch between predefined electrocautery and nerve stimulation waveforms at specified time intervals
  - a. Dynamic impedance tracking → next step
- B. Prototyped circuit from Part. A that can be integrated with an electrosurgical unit and standard cautery probe

## Circuit Development

- Implement pulse width modulation schemes that alternate between cautery and nerve-specific stimulation waveforms, using a high voltage, high frequency ESU as the power source.
- Develop timing circuits to control nerve-specific stimulation pulses.
- Design safe, real-time modulation hardware compatible with existing ESUs.

## Device Prototyping

- Design and test functional PCB layouts
- Integrate hardware with microcontrollers, or if interested, FPGA-based systems

# Working with NeuroGuard for ECE 445

## Opportunities to explore:

- Advanced circuit design and prototyping (power electronics)
- Medical devices
- Surgical/intraoperative considerations (patient safety, energy delivery, neurophysiology, etc.)

## Forward-thinking Team:

- 1st place, 2024 CIMED Health Makeathon
- Finalists, 2025 Global Consortium Healthcare Innovation Challenge
- Semifinalists, 2025 Grainger Engineering Tech Challenge at Chicago Week
- Research and Innovation Pathway Grant, Carle Illinois College of Medicine

**Real-World Impact:** NeuroGuard has the potential to lower global nerve injury rates in mastectomies—especially in resource-limited settings where access to advanced equipment is limited. You'll work closely with medical students, surgeons, and engineering faculty, gaining exposure to translational innovation at the intersection of healthcare and engineering.



# NeuroGuard

*Cutting Edge Precision without Cutting Nerves*

*ECE 445 Senior Design Course Launch, 8/26/25*

